


LC series


H2 GENERATOR DESK

The LC series generators use an electrolytic cell with polymeric membrane (PEM) to produce pure hydrogen. These series features a desiccant cartridge that needs to be replaced or refilled when saturated. The exclusive, electronically-controlled gas/liquid separator, automatic checking for internal leaks whenever starting the unit, and constant control of operating parameters, guarantee maximum safety. Up to 20 units can be connected in parallel. The touch-screen LCD interface provides simple and user-friendly management of all functions on the unit.



- 1 Touch-screen LCD 128x64 pixel
- 2 START/STOP button
- 3 Hydrogen Outlet
- 4 Hydrogen purge
- 5a Water connector for filter
- 5b Water connector for filter and for emptying the tank
- 6 Hydrogen vent
- 7 Oxygen vent
- 8 Desiccant cartridge compartment
- 9 Power connection and switch
- 10 I/O connector: RS232 – RS485 (Optional)
- 11 Cooling fan air intake


Main Applications

- ICP-MS collision gas
- Flame ionization detector feed gas (FID)
- Refilling metal hydride tanks for use with fuel cells

Main Features

- **Available Flow-rates:**
up to 260 cc/min
- **Outlet pressure:**
up to 7 bars
- **Hydrogen purity:**
>99.9999%
- **Drying system:**
Refillable molecular sieve column
- **Internal water tank:**
1.2 litres, with electronic level control
- **Dimensions:**
24x33x31(H)
- **Weight:**
6.5 kg
- **Certification:**
CE, ISO9001

Principle diagram



Hydrogen is produced from distilled water using a polymeric membrane (PEM). No acid or alkaline solutions are used.

The LC series features a desiccant cartridge that needs to be replaced or refilled when saturated.

| Models | LC.H2.100 | LC.H2.140 | LC.H2.180 | LC.H2.260 |
|-------------------------------|-------------------------------|-----------|-----------|-----------|
| General data | | | | |
| Electrolytic cell | PEM technology | | | |
| Drying system | External column dryer | | | |
| H2 purity | >99.9999% ¹ | | | |
| Outlet pressure | 7 bars (100 psi) | | | |
| H2 flow-rate cc/min (Max) | 100 | 140 | 180 | 260 |
| Dimensions | 24x33x31(H) cm | | | |
| Net weight (no water in tank) | 6.5 kg | | | |
| Communication | | | | |
| RS232 | X | X | X | X |
| RS485 | Optional | | | |
| Software functions | | | | |
| Parallel mode | Optional | | | |
| Automatic tank refill | Not available | | | |
| Water | | | | |
| Quality | Deionized, ASTM II, <0.1uS | | | |
| Internal tank capacity | 1.2 l | | | |
| Electrical data | | | | |
| Type of connection | IEC 320-C13 (DESKTOP ADAPTER) | | | |
| Power supply voltage | 100-240 Vac 50/60 Hz | | | |
| Installed power (Max) | 120W | | 160W | |
| Connections | | | | |
| Hydrogen outlet | 1/8" compression fitting | | | |
| Water | Quick release push-in fitting | | | |

¹ Referred to O₂ dew point < -55°C